
# Pharmacoll® & Pharmacoll® *comfort*

## Pharmacoll

Pharmacoll comfort

STERILE	R
---------	---

## Pharmacoll®

**Description:**

**Pharmacol<sup>®</sup>** is a gelatin free hydrocolloid adhesive wound dressing laminated to a semi permeable, bacteria proof, virus proof and waterproof polyurethane film backing individually packed and sterilized by radiation. **Pharmacol<sup>®</sup> comfort** with adhesive polyurethane film border.

**Pharmacol<sup>®</sup> basic** is laminated to PU foam backing providing extra cushion. **Pharmacol<sup>®</sup> thin** has lower absorption capacity, so it could be used during epithelialization phase. **Pharmacol<sup>®</sup>** forms a gel upon absorbing wound exudates providing a moist wound environment, supporting granulation and epithelialization. Gel formation prevents **Pharmacol<sup>®</sup>** from adhering to wounds, and dressing can be removed without causing pain and without damaging the granulation and epithelial tissue.

### Benefits:

- Free of gelatin, minimizing associated allergy and odour.
- Semi-permeable, waterproof, bacteria proof and virus proof polyurethane backing, protects wound from moisture, bacteria and virus.
- High absorption capacity of **Pharmacoll®** binds away excess exudates prolonging dressing change intervals.
- **Pharmacoll®** accelerates healing process by providing ideal environment for formation of granulation tissue.
- Slightly transparent **Pharmacoll®** enables wound inspection.

**Indications:**

- Suitable for moist wound care
- Suitable for light to moderately exudating wounds second degree burns, leg ulcers, pressure ulcers, chronic slow healing wounds and abrasions.
- Used also for blisters treatment.

### Contraindications:

- **Pharmacol<sup>®</sup>** should not be used on third degree burns nor on exposed muscle, bones and tendons nor on infected wounds .
- Discoloration and bulging of the dressing indicate it should be changed.

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