FARMAC-ZABBAN S.p.A.

Indirizzo
Tel. / Fax
e-mail

Via Persicetana, 26 - Calderara di Reno (BO)
+39 051 – 318 411 / +39 051 – 318 472 - 318 447
m.gruppioni@farmaczabban.it company@farmaczabban.it

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LOTTO N. 2

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
Product specification

1. Product name SOL-M™ Blunt Fill Needle

2. Description The SOL-M™ Blunt Fill Needle is designed to replace hypodermic needles currently for withdrawal of medication from rubber- stoppered vials or injection into I.V. Systems and pre-slit septum covering injections sites.

3. Characteristics The SOL-M™ Blunt Fill Needle is used in conjunction with a syringe to penetrate rubber-stoppered medicine vials and pre-slit septum covering injection sites. The cannula is pre-lubricated to reduce septum insertion forces.

4. Intended use The SOL-M™ Blunt Fill Needle is used in conjunction with a syringe as an additive device for aspiration from multi-dose medicine vials or injection into I.V. Systems and pre-slit septum covering injection sites.

5. Instructions for use See current IFU


6. Sizes and REF numbers

REF	Size
110021	18GX1"
110022	18GX1 1/2"
110023	18GX2"

Technical information

		Component name	Material
1. List of Materials		Needle hub	PP: H1500
		Needle tube	Stainless steel: SUS304
		Needle cap	PP: RP344RK
		Adhesive	Epoxy: YD-128
		Needle Lubricant	Silicon oil: Xilikang 201-350cst
2. Latex free		YES	
3. PHT / DEHP / PVC / BPA free		YES	
4. Shelf life		5 years	
5. Sterilization method		Sterilized using Ethylene Oxide	
6. Packaging specification	6.1 Sales unit	100	Units per box
		1000	Units per case

7. Technical Drawing


- 1. Needle hub
- 2. Needle tube
- 3. Needle cap

Quality and Regulatory information



- | | |
|---------------------------|--|
| 1. Quality certificate | Quality Management System according ISO 13485 |
| 2. Product classification | Class IIa according to Annex IX of MDD 93/42/EEC |

The product is compliant with the following standards and regulations:

	Document reference	Title
3. List of standards	ISO 7864:2016	Sterile hypodermic needles for single use -- Requirements and test methods
	ISO 9626:2016	Stainless steel needle tubing for the manufacture of medical devices -- Requirements and test methods
	ISO 80369-7:2016	Small-bore connectors for liquids and gases in healthcare applications -- Part 7: Connectors for intravascular or hypodermic applications
	ISO10993-4:2002/Amd 1:2006	Biological evaluation of medical devices -- Part 4: Selection of tests for interactions with blood
	ISO10993-5:2009	Biological evaluation of medical devices -- Part 5: Tests for in vitro cytotoxicity
	ISO10993-7:2008/Cor 1:2009	Biological evaluation of medical devices -- Part 7: Ethylene oxide sterilization residuals
	ISO10993-10:2010	Biological evaluation of medical devices -- Part 10: Tests for irritation and skin sensitization
	ISO10993-11:2006	Biological evaluation of medical devices -- Part 11: Tests for systemic toxicity
	ISO 15223-1:2016	Medical devices -- Symbols to be used with medical device labels, labelling and information to be supplied -- Part 1: General requirements
	ISO 11607-1:2006/Amd 1:2014	Packaging for terminally sterilized medical devices -- Part 1: Requirements for materials, sterile barrier systems and packaging systems
	ISO 11607-2:2006/Amd 1:2014	Packaging for terminally sterilized medical devices -- Part 2: Validation requirements for forming, sealing and assembly processes
	EN1041:2008+A1:2013	Information supplied by the manufacturer of medical devices

Non-pyrogenic/Non-toxic
Pyrogenfrei/Nicht toxisch
Non pirogénico/No tóxico
Non pyrogène/Non toxique
Apyrogeno/Non tossico
Não pirogénico/Não tóxico
Pyrogenfri/Ej toxisk
Pyrogenfri/Myrkytön
Pyrogenfri/lkke-toksisk
Pyrogenfri/lkke-toksisk
Pyrogenfri/Niet-toxisch
Níepirogennj/Nietolszysznen
Mij puropetogyno/Mij toxičko

Blunt Fill Needle w/5 Micron Filter

The Sol-M™ Blunt Fill Needle with 5 Micron Filter is used for drawing up medication from an ampule. This should be used in combination with any standard Luer Lock Syringe while drawing up medication from an ampule.

1. Remove the Blunt Fill Needle from its package and connect to the syringe.(See Figure 1 & 2)
2. Remove the needle cover.
3. Blunt Fill Needle with Filter: Insert the needle into the center of the vial stopper at a 90° angle. This will help reduce the risk of rubber fragments from the side of the vial stopper contaminating the medication. Draw up medication from inverted glass ampoule (without rubber stopper). (See Figure 3)
4. Once medication has been prepared and is in the syringe ready for administration, remove the Blunt Fill Needle with Filter and dispose of in a sharps container.(See Figure 4)

The Blunt Fill Needle is for single use only. Reuse will dull the needle, remove the silicone coating and increase the chance of coring and rubber fragmentation.

Caution:

This device is sterile unless the package has been damaged or opened.
This device is a single use item only.

This device is a single use item only.
Sol-M™ Blunt Fill Needle with 5 Micron Filter is to be used for draw up of medication from an ampule only, not for patient injection.

Federal Law (USA) restricts this device to sale by or on the order of a physician. Follow standard infection control/precautions as a specified by the Centers Disease Control and Prevention (USA) or local equivalent.

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

Blunt Fill Needle

The Sol-M™ Blunt Fill Needle should be used in combination with any standard Luer Lock Syringe for draw up of medication.

1. Remove the Blunt Fill Needle from its package and connect to the syringe. (See Figure 1 & 2)
2. Remove the needle cover.
3. Blunt Fill Needle: Insert the needle into the center of the vial stopper at a 90° angle. This will help reduce the risk of rubber fragments from the side of the vial stopper contaminating the medication. (See Figure 3)
4. Once medication has been prepared and is in the syringe ready for administration, remove the Blunt Fill Needle dispose of in a sharps container. (See Figure 4)

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Caution:

This device is sterile unless the package has been damaged or opened. This device is a single use item only.
Sol-M™ Blunt Fill Needle is designed for penetration into vials and accessing IV Bags, not for patient injection.

Federal Law (USA) restricts this device to sale by or on the order of a physician. Follow standard infection control precautions as a specified by the Centers Disease Control and Prevention (USA) or local equivalent.

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TO:	To whom it may concern
FROM:	Ligang Xue, Director, Product development
DATE:	2016-05-11

Sol-Millennium Medical Inc. is the legal manufacturer for all SOL-CARE™, SOL-M™ and SOL-GUARD™ branded products. We hereby declare, that they meet provisions of the Council Directive 93/42/EEC, the products are Phthalate free, as well as are packed in a pyrogen-free, latex-free (natural rubber-free) and recyclable packaging.

All products in compliance are listed in the following table:

1. SOL-CARE™ Safety Needle;
2. SOL-CARE™ Luer Lock Syringe with Safety Needle;
3. SOL-CARE™ Insulin Safety Syringe U-100 Insulin Only;
4. SOL-CARE™ TB Safety Syringe with Fixed Needle;
5. SOL-CARE™ Luer Lock Safety Syringe with Exchangeable Needle;
6. SOL-CARE™ Luer Lock Safety Syringe without Needle;
7. SOL-GUARD™ TB Safety Syringe with fixed needle;
8. SOL-GUARD™ Insulin Safety Syringe U-100 Insulin Only;
9. SOL-M™ Blunt Fill Needle W/5 Micron Filter;
10. SOL-M™ Blunt Fill Needle;
11. SOL-M™ Oral Dispensing Syringe Amber with Tip Cap;
12. SOL-M™ Oral Dispensing Syringe Purple Plunger with Tip Cap;
13. SOL-M™ Oral Dispensing Syringe Clear with Tip Cap;
14. SOL-M™ Insulin Syringe with Fixed Needle;
15. SOL-M™ Tuberculin Syringe;
16. SOL-M™ Luer Lock Syringe without Needle;
17. SOL-M™ Luer Lock Syringe with Exchangeable Needle;
18. SOL-M™ Hypodermic Needle;
19. SOL-M™ Bottle Adapter
20. SOL-CARE™ Safety Blood Collection Tube Holder Set with Vacuum Activation Tube

Ligang Xue,
Director, Product Development
Sol-Millennium Medical Group

SOL-Millennium Medical Inc.
1735 North Brown Road, Suite 120, Lawrenceville, Georgia, 30043 USA
Office 404.973.2200 Fax 404.973.2212
www.sol-m.com