
Cochlear Implants

MED⁹**EL**

Impianto al tronco encefalico Mi1200 SYNCHRONY

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






AW31686-6.0 (Italian)


hearLIFE

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SYNCHRONY ABI

Typical dimensions in mm
Dimensioni tipiche in mm

Instructions for use

Mi1200 SYNCHRONY Auditory Brainstem Implant

Device description

The Mi1200 SYNCHRONY Auditory Brainstem Implant (hereafter referred to as SYNCHRONY ABI) is the implantable part of the MED-EL ABI System and can only be used together with compatible MED-EL external components. The device consists of a stimulator, a coil with a removable magnet within its centre, a reference electrode, an EAP reference electrode and an active electrode permanently attached to the stimulator. This device is intended to be implanted by adequately trained and experienced surgeons only.

The implant offers a stimulation mode and a telemetry mode. Stimulation sequences of biphasic and triphasic pulses can be delivered sequentially or simultaneously on two or more channels. In telemetry mode the device allows a functional check about the technical status of the implant including communication over the transcutaneous link as well as the assessment of the electrode impedances and recording of the electrically evoked compound action potential.

The implant has a mass of 7.6 g (typical weight).

For principal dimensions of the implant refer to the drawings on the previous page.

The volume of the implant without electrode is 3.7 cm³.

Following materials are in direct contact with human tissue: medical grade silicone, platinum, iridium, polyester and parylene c.

Purpose of the device

The purpose of the device is the perception of environmental sound and potential for improvement of communicational abilities.

Performance characteristics

- Output characteristics of a stimulation signal on a 1 kOhm resistor:
Maximum current amplitude: Median value = 1250 µA, range = 500 µA
Maximum pulse width: Median value = 203.8 µs, range = 8.2 µs
- The impedance measurement accuracy is typically better than 5 %.
- The implant is MR Conditional for scanner field strengths of 0.2T, 1.0T and 1.5T. For more information please refer to the Medical Procedures Manual.
- There are no default factory settings of the implant system.

- Proper functioning of the implantable part of the system can be checked by performing telemetry (refer to the MED-EL application software user manual).

Specification and characteristics for each lead and electrode array

- The implant has 24 independent current sources stimulating 12 independent electrode channels in monopolar mode.
- The electrode is made of medical grade silicone, platinum (electrode contacts), platinum/iridium (90/10) wires and polyester.
- The electrode does not deliver any medicinal substances.
- Physical dimensions of the electrodes:

Electrode type	Length of the electrode lead (typical value including electrode array)	Dimensions of the electrode array		Geometric surface area per channel		Distance	
		Length and width of electrode array *	Thickness of electrode array *	of smallest stimulating electrode contact **	of largest stimulating electrode contact **	between contacts *	between most proximal and most distal contact *
ABI	111.4 mm	5.5 × 3	0.6	0.237	0.237	1	4.5

* typical value, mm

** typical value, mm²

The implant does not have any connector.

Intended use

The SYNCHRONY ABI is used for electrical stimulation of the cochlear nucleus (CN) via an implanted stimulator and a specially designed electrode array to evoke auditory sensations in patients with non-functional cochlear nerves.

Indications

The SYNCHRONY ABI is indicated for patients 15 years of age or older who have been diagnosed with Neurofibromatosis Type 2 (NF2) when both cochlear nerves are non-functional, or are anticipated to be rendered non-functional by presence or removal of tumour. The surgical procedure for the implantation of the device shall be concurrent with tumour removal surgery.

Additionally, the SYNCHRONY ABI is indicated for patients of age 12 months and older who cannot benefit from cochlear implant due to non-functioning auditory nerve. The cause may be congenital due to auditory nerve aplasia or auditory nerve hypoplasia or accrued during the life span due cochlear nerve disruption due to head trauma, non-NF2 tumour or severe cochlear ossification.

To obtain the optimal benefit from the implant, the prospective implant users and their families shall be highly motivated and have realistic expectations about the expected benefit of the implant and shall understand the importance of returning to the implant centre for regular audio processor programming, assessment sessions and training.

Contraindications

A patient must not be implanted,

- if the individual is known to be intolerant of the materials used in the implant (including medical grade silicone, platinum, iridium, polyester and parylene c);
- if the cochlear nerve can be preserved during tumour removal surgery;
- if deafness is attributed to central damage of the upper auditory pathway;
- if anatomical abnormalities or surgical complications which might prevent placement of the SYNCHRONY ABI Active Electrode Array are present;
- if an anatomic situation is present, that would prevent appropriate placement of the stimulator housing in the bone of the skull, the implant shall be carefully considered prior to surgery;
- if the medical team decides that the medical situation of the patient is critical due to multiple tumours, hydrocephalus and other possible adverse status;
- if the psychological status of the patient is unstable or if the patient has unrealistic expectations;
- if the physical condition of the patient prohibits ABI surgery.

Undesirable side effects – Risks related to the implant

Possible postoperative side effects include the following: dizziness, increased vertigo, delay of healing of the scar, impairment of the sense of taste, potential for swallowing difficulties, numbness, increased tinnitus, stimulation of the facial nerve, temporary pain and uncomfortable sounds during stimulation. Though rare, there is also a possibility of postoperative device failure or a decrease in device performance. Some postoperative complications may require revision surgery.

Sterility

The implant has been subjected to a validated ethylene oxide sterilisation process and is supplied in sterile packaging. Once the sterile packaging has been opened, the implant cannot be resterilised. Do not use if sterile packaging is damaged. The implant is for single use only. Do not remove from sterile packaging until required.

Storage, shipment and disposal

The sterilised implant may only be shipped¹ between -20 °C (-4 °F) and +55 °C (+131 °F) and stored inside the implant box at room temperature. Each device must be implanted before the "use by" date specified on the package. Packaging² should be disposed of in accordance with local legislation.

Information about use – General precautions and warnings

- The device must not be altered and must only be used as intended.
- Expected performance with the SYNCHRONY ABI cannot be accurately predicted. Prospective SYNCHRONY ABI users must be properly counselled prior to surgery. The prospective implant users and their families shall be highly motivated and have realistic expectations about the expected benefit of the implant.
- Long-term damage to neural tissue following continuous chronic electrical stimulation has not been observed with auditory brainstem implants.
- Sterility of the implant must be ensured at all times.
- The implant must never be dropped onto a hard surface or be held only by the electrode; damage to the implant or electrodes during the operation will invalidate the warranty.
- Device failure may occur due to mechanical damage of the implanted parts, which could result from a blow to the head, or due to electronic or other technical failure of the implant. Reimplantation of the device is required in these cases.
- Make sure that the length of the electrode lead allows placement of the electrode paddle without compressing or extending the lead.
- The platinum bridge on the overmould at implant housing must be removed before closing the wound.
- The SYNCHRONY ABI Active Electrode can be prefixed onto the brainstem surface with a piece of hemostyptic material and fibrin glue to prevent postsurgical displacement. The main fixation is achieved by gluing the electrode lead to the rostral surface of the cerebellum. Useful for this procedure is collagen foam covered with instant glue (e.g. fibrin glue).
- If necessary, the looping of the electrode lead shall be performed outside of the dura within the mastoid cavity.
- After implantation, the dura shall be closed in a watertight fashion. The entry of the electrode lead into the dura is sealed with fibrin glue or collagen foam covered with instant glue (e.g. fibrin glue).
- Particular attention should be paid to paediatric patients with developmental challenges (e.g. Cornelia-de-Lange syndrome), as they may present an increased risk of suffocation and/or laryngeal injury from attempting to swallow the coil and/or other external parts.
- Most water sports should not cause any problem as long as the external parts of the implant system are removed. If headgear or face mask are worn, care must be taken to ensure that the strap is not too tight over the site of the implant. In case the MED-EL implant user wants to dive, the user should consult an experienced physician about the possibilities and personal restrictions when performing water sports, especially in the case of SCUBA diving. The implant is robust against pressure changes which occur during SCUBA diving to depths up to 50m (165 ft).

Surgical precautions and warnings – Risks related to surgery

- If a tumour removal surgery is necessary or planned, implantation of the device should be performed during the same surgical intervention. Additional time required for implantation of the device increases the standard risk associated with tumour removal surgery.
- If previous irradiation of the tumour has taken place, the patient's evaluation should include possible structural deformation caused by such irradiation. Large structural deformation may prevent proper SYNCHRONY ABI placement and should be evaluated by the surgeon.
- The normal risks of surgery and general anaesthesia are applicable. Primary surgical risks include the following: infection, inflammation, swelling, necrosis, haematoma, leakage of CSF, damage to the facial nerve, pain, scarring of the wound, skin irritation, swallowing difficulties and complications related to general anaesthesia. Additionally, meningitis³ can be a rare postoperative complication, but has the potential to be serious. The risk of meningitis may be reduced, for example by vaccination, antibiotic use, and surgical technique.
- Intensive intraoperative monitoring is necessary during acoustic neuroma surgery and the implantation of the SYNCHRONY ABI to control and protect nerve structures in and around the brainstem. Cochleo-vestibular (VIII.), facial (VII.) and caudal cranial (IX., X.) nerve structures should be monitored during surgery.
- The anaesthetist shall be familiar with the monitoring, detection and management of air embolism and should also be familiar with anaesthetics connected with electrophysiological measurements.
- Prior to the implantation of the SYNCHRONY ABI, EABR measurements should be carried out.
- Prophylactic use of antibiotics is recommended for all patients unless medically contraindicated.
- In order to achieve good magnetic holding power and optimal coupling the distance between the lateral side of the implant and the surface of the skin (with hair) shall not exceed 6 mm.
- The serial number of the implant must be visible on the implant before fixing it in place.
- The implant must be immobilised in a flat stimulator bed drilled in the temporal bone. The electrode lead should be placed in a ramp-like bony channel without sharp edges to protect it against postoperative movement and excessive mechanical impact. The anterior stimulator edge should not be recessed to a depth more than 2 mm.
- Additional immobilisation of the implant needs to be done (e.g. with sutures). It should be done in such a way that there will be no postoperative movement. Continuous movement may result in mechanical fatigue and subsequent premature failure of electrical connections.

- The electrode array should be inserted into the lateral recess (LR) of the IV. ventricle without squeezing either the electrode array or the electrode lead and without touching the electrode contacts.
- Do not place the sutures directly over the electrode lead.
- Inaccurate placement of the electrode array may impair acoustic perceptions with the device and may necessitate additional surgery. Improper fixation or placement of the electrode array may also result in premature failure of the implant. Localised non-auditory sensation such as tingling in the head and neck, upper extremities, torso, and lower extremities may take place. Non-auditory sensation could be localised in the ipsilateral or contralateral side of the implant.
- Improper fixation or placement of the SYNCHRONY ABI Electrode Array may also result in dislocation (e.g. in case of large lateral recess) of the SYNCHRONY ABI Electrode Array or premature failure of the implant. Revision surgery may be required in these cases, if possible.
- Device failure may also occur due to mechanical trauma to the implanted parts, such as a blow to the head, or due to electronic or other technical failures of the implant. Reimplantation may be required in these cases, if possible.
- Good physical and thus stable electrical contact between stimulation reference electrode and surrounding tissue is essential for electrical stimulation. Therefore do not place any fixation sutures directly over the reference electrode and do not recess the stimulator too deeply to avoid any air gap over the reference electrode.
- Only MED-EL approved and recommended surgical instruments must be used during the insertion process.
- The implant contains a strong magnet. Never use magnetic surgical tools.
- The initial fitting of the audio processor is done typically six to eight weeks after surgery. The fitting shall be performed in the intensive care unit of the hospital or in the presence of an intensive care crash cart.
- For patients that are not able to give appropriate feedback (e.g. paediatric patients) it is recommended to perform EABR measurement under general anaesthesia in an Intensive Care Unit (ICU) before initial fitting. During the EABR recordings every electrode needs to be checked for any non-auditory perception.
- In case of a magnet exchange surgery MED-EL recommends the use of the following tools: The Magnet Removal Tool Ms050206 or the Magnet Insertion Tool Ms050205 can be used in combination with the Non-Magnetic Spacer Ms010107 or the Replacement Magnet Ms010108. Please refer to the applicable Instructions for Use of these devices for further surgical precautions and warnings.
- After removal of the Non-Magnetic Spacer Ms010107, make sure that the Replacement Magnet Ms010108 is inserted to re-establish full functionality of the implant.
- Other risks after surgery may be avoided by following the instructions in the applicable MED-EL audio processor user manual and the MED-EL application software user manual.

Interference with other equipment, robustness of the device in special medical or diagnostic environments

For all other medical procedure-related safety recommendations and guidelines, including MRI, please refer to the Medical Procedures Manual.

Explantation

The implant may become non-functional, either by accident or due to medical or technical reasons. In this case it is strongly recommended to explant the device.

If for any reason the device is not used anymore, it is strongly recommended to explant the device. If an explantation is not performed functional checks of the implant on a regular basis are strongly recommended.

If possible, the device should be removed without damaging or cutting it.

MED-EL recommends conducting a functionality check (telemetry measurement) as well as taking an x-ray image of the device in implanted condition prior to the explantation of the device.

Further details regarding the device condition as found during explantation surgery are requested within MED-EL's "Device Explant Report Form" (AW36523), which is contained in MED-EL's "Explant Kit" (PN04175).

Please make sure that the "Device Explant Report Form" is completed and provided along with functionality check results, x-ray images and the explanted device.

Returning explanted devices

- After the device has been surgically removed, follow the cleaning and disinfection procedures established at the explantation site, avoiding damage to the implant if possible. Follow locally established procedures for potentially bio-hazardous material at all times.
- The device is to be returned to MED-EL in the Returned Implant Kit. Follow the enclosed packaging instruction.
- The device is returned to:
MED-EL Elektromedizinische Geräte GmbH
Worldwide Headquarters
Safety adviser for the transport of dangerous goods
Fürstenweg 77a
6020 Innsbruck
Austria

Warranty

Please refer to the accompanying Warranty Statement for information on our warranty provisions.

Symbols

CE mark. First applied in 2014

Single-use device. Do not reuse!

Catalogue number

Serial number

Caution! Consult accompanying documents

Sterilised using Ethylene Oxide

"Use by" date

Manufacturing date

Manufacturer

MR Conditional

Help and assistance are always available from your local office.
Please refer to the accompanying Contact Sheet for your local office.
Please visit us at www.medel.com

-
- 1 For shipping, the Implant Box shall be packed into a protective padded cardboard box (or similar).
 - 2 The cardboard and plastic implant packaging (PETG and Tyvek) are manufactured from recyclable materials.
 - 3 A paper reporting on the pathomechanisms, clinical symptoms, conservative and surgical treatments in cases of meningitis, published by Arnold et al (ORL 2002;64:382-389), may be useful additional reading.

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