Promus ELITE™

MONORAIL™

Everolimus-Eluting Platinum Chromium
Coronary Stent System

Directions for Use	2
Instrucciones de uso	11
Mode d'emploi	21
Gebrauchsanweisung	31
Istruzioni per l'uso	41
Gebruiksaanwijzing	51
Instruções de Utilização	61

TABLE OF CONTENTS

1.0 WARNING	3	7.0 DRUG INFORMATION	7
2.0 DEVICE DESCRIPTION	3	7.1 Mechanism of Action	7
Table 2.1 Promus ELITE™ Stent System Product Description	3	7.2 Drug Interactions	7
2.1 User Information	3	7.3 Carcinogenicity, Genotoxicity, and Reproductive Toxicology	7
2.2 Device Component Description	3	8.0 ADVERSE EVENTS	7
Contents	3	9.0 HOW SUPPLIED	8
2.3 Drug Eluting Coating Description	3	Handling and Storage	8
2.3.1 Everolimus	3	10.0 OPERATIONAL INSTRUCTIONS	8
Figure 2.1. The Chemical Structure of Everolimus	3	10.1 Inspection Prior to Use	8
2.3.2 Primer Polymer and Drug Matrix Copolymer Carrier	3	10.2 Materials Required (not included in Stent Delivery System package)	8
Figure 2.2. The Chemical Structure of PBMA	4	10.3 Preparation	8
Figure 2.3. The Chemical Structure of PVDF-HFP	4	10.3.1 Packaging Removal	8
Table 2.2. Promus ELITE Stent System Product Matrix and Everolimus Content	4	10.3.2 Guidewire Lumen Flush	9
3.0 INTENDED USE/INDICATIONS FOR USE	4	10.3.3 Balloon Preparation	9
4.0 CONTRAINDICATIONS	5	10.3.4 Delivery Procedure	9
5.0 WARNINGS	5	10.3.5 Deployment Procedure	9
6.0 PRECAUTIONS	5	10.3.6 Removal Procedure	9
6.1 General Precautions	5	10.4 Post-Deployment Dilatation of Stented Segments	10
6.2 Stent System Handling (also see Section 10.0, Operational Instructions)	5	11.0 IN VITRO INFORMATION	10
6.3 Stent Placement	5	11.1 Promus ELITE Stent System Compliance	10
6.3.1 Preparation	5	Table 11.1. Typical Promus ELITE Stent System Compliance	10
6.3.2 Placement	5	Table 11.2. Typical Promus ELITE Stent System Compliance	10
6.4 Stent System Removal – Pre-deployment	5	12.0 CLINICAL TRIALS	10
6.5 Stent System Removal – Post-deployment	6	WARRANTY	10
6.6 Post-Procedure	6		
6.7 Brachytherapy	6		
6.8 Magnetic Resonance Imaging (MRI)	6		
3.0 Tesla Temperature Information	6		
1.5 Tesla Temperature Information	6		
Image Artifact Information	6		
6.9 Individualization of Patient Treatment	6		
6.10 Drug Interactions	7		
6.11 Pregnancy	7		
6.12 Use of Multiple Stents	7		

Promus ELITE™

MONORAIL™

Everolimus-Eluting Platinum Chromium Coronary Stent System

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

1.0 WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. **DO NOT REUSE, REPROCESS OR RESTERILIZE.** Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

STERILE — DO NOT RESTERILIZE — SINGLE USE ONLY

2.0 DEVICE DESCRIPTION

The Promus ELITE Everolimus-Eluting Coronary Stent System (Promus ELITE Stent System) is a device/drug combination product composed of two components: a device (coronary stent system) and a drug product (a formulation of everolimus contained in a polymer coating). The characteristics of the Promus ELITE Stent System are described in Table 2.1 Promus ELITE Stent System Product Description.

Table 2.1 Promus ELITE Stent System Product Description

	Promus ELITE Monorail Stent Delivery System
Available Stent Lengths (mm)	8, 12, 16, 20, 24, 28, 32, 38*
Available Stent Diameters (mm)	2.25*, 2.50, 2.75, 3.00, 3.50, 4.00
Stent Material	Platinum Chromium Alloy (PtCr)
Drug Product	A conformal coating of a polymer carrier loaded with 100 µg/cm² everolimus applied to the stent with a maximum nominal drug content of 243.0 µg on the largest stent (4.00 x 38 mm).
Delivery System Effective Length	144 cm
Delivery System Ports	Single access port to inflation lumen. Guidewire exit port is located approximately 26 cm from tip. Designed for guidewire ≤0.014 inches (0.36 mm)
Average stent length change at nominal diameter	2.25 – 4.00 mm; 0.1 - 1.5 mm
Stent Delivery Balloon	A balloon, with two radiopaque balloon markers, nominally placed 0.4 mm (0.016 inches) beyond the stent at each end.
Balloon Inflation Pressure	Nominal Inflation Pressure: 11 atm – 1117 kPa Rated Burst Pressure: 18 atm – 1827 kPa for stent diameters 2.25 – 2.75 (mm) and 16 atm – 1620 kPa for stent diameters 3.00 – 4.00 (mm) diameters

	Promus ELITE Stent Delivery System
Guide Catheter Inner Diameter	≥0.056 inches (1.42 mm)
Catheter Shaft Outer Diameter	2.1F (0.70 mm) proximal and 2.7F (≤0.95 mm) distal.
Stent Strut Thickness (including coating)	2.25 – 3.50 mm: 0.093 mm, 4.00 mm: 0.098 mm

*38 mm length is not available in 2.25 mm diameter sizes.

2.1 User Information

Only Physicians who have received adequate training should perform implantation of the stent.

2.2 Device Component Description

The Promus ELITE Stent component is made from a platinum chromium alloy mounted onto a Monorail Delivery System. The Promus ELITE Stent is available in 4 stent models each designed for specific diameters as follows:

- Small Vessel (SV): 2.25 mm
- Small Workhorse (SWH): 2.50, 2.75 mm
- Workhorse (WH): 3.00, 3.50 mm
- Large Vessel (LV): 4.00 mm

Contents

Qty	Material
One (1)	Promus ELITE Monorail Stent System
One (1)	Flushing needle with luer fitting

2.3 Drug Eluting Coating Description

The Promus ELITE Stent is a stent with a drug/polymer coating. The coating is comprised of two layers: the inner layer consists of a polymer which is a primer for improving the adhesion of the outer layer; the outer layer is a polymer matrix that contains an active pharmaceutical ingredient.

See Sections 2.3.1 and 2.3.2 for descriptions of the drug and polymers, respectively.

2.3.1 Everolimus

The active pharmaceutical ingredient in the Promus ELITE Stent is everolimus. The everolimus chemical name is 40-O-(2-hydroxyethyl)-rapamycin and its chemical structure is provided in Figure 2.1. This is the same active pharmaceutical ingredient as is used in PROMUS™ (Xience V™), PROMUS Element™, PROMUS Element™ Plus and Promus Premier™.


Figure 2.1. The Chemical Structure of Everolimus

2.3.2 Primer Polymer and Drug Matrix Copolymer Carrier

The Promus ELITE Stent contains a primer polymer layer, PBMA, poly (n-butyl methacrylate), that functions as an adhesion promoter between the bare metal and the drug matrix layer. The chemical structure of PBMA is provided in Figure 2.2.


Figure 2.2. The Chemical Structure of PBMA

The drug matrix layer contains a semi-crystalline random copolymer, PVDF-HFP, poly (vinylidene fluoride-co-hexafluoropropylene), blended with everolimus. The chemical structure of PVDF-HFP is provided in Figure 2.3.


Figure 2.3. The Chemical Structure of PVDF-HFP

Table 2.2. Promus ELITE™ Stent System Product Matrix and Everolimus Content

Product Code MR	Nominal Expanded Stent Inner Diameter (mm)	Nominal Unexpanded Stent Length (mm)	Nominal Everolimus Content (µg)
H7493941308220	2.25	8	38.2
H7493941308250	2.50	8	39.3
H7493941308270	2.75	8	39.3
H7493941308300	3.00	8	42.6
H7493941308350	3.50	8	42.6
H7493941308400	4.00	8	57.3
H7493941312220	2.25	12	57.3
H7493941312250	2.50	12	61.1
H7493941312270	2.75	12	61.1
H7493941312300	3.00	12	60.7
H7493941312350	3.50	12	60.7
H7493941312400	4.00	12	81.5
H7493941316220	2.25	16	72.7
H7493941316250	2.50	16	78.5
H7493941316270	2.75	16	78.5
H7493941316300	3.00	16	84.8
H7493941316350	3.50	16	84.8
H7493941316400	4.00	16	105.7
H7493941320220	2.25	20	91.8
H7493941320250	2.50	20	95.8
H7493941320270	2.75	20	95.8

Product Code MR	Nominal Expanded Stent Inner Diameter (mm)	Nominal Unexpanded Stent Length (mm)	Nominal Everolimus Content (µg)
H7493941320300	3.00	20	102.9
H7493941320350	3.50	20	102.9
H7493941320400	4.00	20	129.9
H7493941324220	2.25	24	107.2
H7493941324250	2.50	24	113.2
H7493941324270	2.75	24	113.2
H7493941324300	3.00	24	121.1
H7493941324350	3.50	24	121.1
H7493941324400	4.00	24	154.1
H7493941328220	2.25	28	126.3
H7493941328250	2.50	28	130.6
H7493941328270	2.75	28	130.6
H7493941328300	3.00	28	139.2
H7493941328350	3.50	28	139.2
H7493941328400	4.00	28	178.4
H7493941332220	2.25	32	145.5
H7493941332250	2.50	32	152.3
H7493941332270	2.75	32	152.3
H7493941332300	3.00	32	163.3
H7493941332350	3.50	32	163.3
H7493941332400	4.00	32	202.6
H7493941338250	2.50	38	178.4
H7493941338270	2.75	38	178.4
H7493941338300	3.00	38	193.5
H7493941338350	3.50	38	193.5
H7493941338400	4.00	38	243.0

3.0 INTENDED USE/INDICATIONS FOR USE

The Promus ELITE Everolimus-Eluting Platinum Chromium Coronary Stent System is indicated for improving coronary luminal diameter in patients with symptomatic ischemic heart disease, including patients with acute myocardial infarction and patients with concomitant diabetes mellitus due to discrete de novo native coronary artery lesions. The Promus ELITE Everolimus-Eluting Platinum Chromium Coronary Stent System is also indicated for treatment of patients presenting with:

- Coronary bifurcation lesions
- Coronary artery ostial lesions
- Unprotected left main coronary artery lesions
- Coronary artery total occlusion lesions
- In-stent restenosis in coronary artery lesions

The treated lesion length should be less than the nominal stent length (8 mm, 12 mm, 16 mm, 20 mm, 24 mm, 28 mm, 32 mm and 38 mm) with a reference vessel diameter of 2.25 mm - 4.00 mm.

4.0 CONTRAINDICATIONS

Use of the Promus ELITE™ Stent System is contraindicated in patients with the following:

- Known hypersensitivity to platinum, the platinum chromium alloy (or similar alloy types such as stainless steel).
- Known hypersensitivity or contraindication to everolimus or structurally related compounds.
- Known hypersensitivity to the polymer or its individual components (see details in Section 2.2.2, Primer polymer and drug matrix copolymer carrier).
- Known severe reaction to contrast agents that cannot be adequately pre-medicated prior to the Promus ELITE Stent placement procedure.

Coronary artery stenting is contraindicated for use in the following:

- Patients who cannot receive recommended anti-platelet and/or anticoagulant therapy.
- Patients judged to have a lesion that prevents complete inflation of an angioplasty balloon or proper placement of the stent or delivery device.

5.0 WARNINGS

- This product should not be used in patients who are not likely to comply with recommended anti-platelet therapy.
- To maintain sterility, the inner packages should not be opened or damaged prior to use.
- The use of this product carries the risks associated with coronary artery stenting, including stent thrombosis, vascular complications, and/or bleeding events.
- Patients with known hypersensitivity to stainless steel, platinum, chromium, iron, nickel or molybdenum, may suffer an allergic reaction to this implant.

6.0 PRECAUTIONS

6.1 General Precautions

- Only physicians who have received adequate training should perform implantation of the stent.
- Stent placement should only be performed at hospitals where emergency coronary artery bypass graft (CABG) surgery is readily available.
- The Heart Team, potentially made up of clinical or non-invasive cardiologists, cardiac surgeons and interventional cardiologists, can be utilized to provide a balanced, multidisciplinary decision-making process on optimal patient care for complex pathologies in accordance with current ESC and/or other local guidelines.
- Subsequent stent blockage may require repeat dilatation of the arterial segment containing the stent. The long-term outcome following repeat dilatation of previously implanted stents is not well characterised.
- Consideration should be given to the risks and benefits of use in patients with a history of severe reaction to contrast agents.
- Do not expose the stent delivery system to organic solvents such as alcohol or detergents.
- Care should be taken to control the position of the guide catheter tip during stent delivery, deployment, and balloon withdrawal. Ensure balloon is fully deflated before delivery system withdrawal. Larger and longer balloons will take more time to deflate than smaller and shorter balloons. Deflation time is ≤30 seconds. Before withdrawing the stent delivery system visually confirm complete balloon deflation under fluoroscopy. Failure to do so may cause increased SDS withdrawal forces and result in guide catheter movement into the vessel and subsequent arterial damage.
- Orally administered everolimus combined with cyclosporine is associated with increased serum cholesterol and triglyceride levels.

6.2 Stent System Handling (also see Section 10.0, Operational Instructions)

- For single use only. Do not resterilize or reuse this product. Note product "Use By" date. (see Section 1.0, Warning)
- The premounted Promus ELITE Stent and its delivery system are designed for use as a unit. The stent is not to be removed from its delivery balloon. The stent is not designed to be crimped onto another balloon. Removing the stent from its delivery balloon may damage the stent and coating and/or lead to stent embolization.

- Prior to angioplasty, carefully examine all equipment to be used during the procedure including the dilatation catheter to verify proper function.
- Special care must be taken not to handle or in any way disrupt the stent position on the delivery balloon. This is most important during catheter removal from packaging, placement over guidewire, and advancement through hemostasis valve adapter and guide catheter hub.
- Excessive manipulation or handling may cause coating damage, contamination, or dislodgment of the stent from the delivery balloon.
- Use only the appropriate balloon inflation media (see Section 10.3.3, Balloon Preparation). Do not use air or any gas medium to inflate the balloon.
- In the event the Promus ELITE Stent is not deployed, do not reuse the product and contact your local Boston Scientific Representative for return information.

6.3 Stent Placement

6.3.1 Preparation

- Do not prepare or pre-inflate balloon prior to stent deployment other than as directed. Use the balloon purging technique described in Section 10.3.3, Balloon Preparation.
- If unusual resistance is felt at any time during lesion access before stent implantation, see **Stent System Removal - Pre-deployment** for directions.
- An unexpanded stent should be introduced into the coronary arteries one time only. An unexpanded stent should not be used subsequent to being moved in and out through the distal end of the guide catheter as stent or coating damage or stent dislodgment from the balloon may occur.

6.3.2 Placement

- The vessel should be pre-dilated with an appropriately sized balloon. Failure to do so may increase the risk of placement difficulty and procedural complications.
- Do not expand the stent if it is not properly positioned in the vessel (see Precautions, Section 6.4, Stent System Removal - Pre-deployment).
- Balloon pressures should be monitored during inflation. Do not exceed the rated burst pressure as indicated on the product label (see Tables 11.1 and 11.2, Typical Promus ELITE Stent System Compliance). Use of pressures higher than specified on the product label may result in a ruptured balloon or shaft. This may result in potential intimal damage, dissection or vessel rupture.
- The stent inner diameter should approximate 1.1 times the distal reference vessel diameter.
- Placement of the stent has the potential to compromise side branch patency if stenting near a side branch.
- Implanting a stent may lead to dissection of the vessel distal and/or proximal to the stented portion and may cause acute closure of the vessel requiring additional intervention (e.g., CABG, further dilation, placement of additional stents, or other).
- When treating multiple lesions, the distal lesion should generally be stented first, followed by stenting of the more proximal lesion(s). Stenting in this order avoids the requirement to cross the proximal stent when placing the distal stent and reduces the chances of the stent dislodgment.
- When treating coronary bifurcation lesions, care must be exercised to access the secondary vessel via the repeating geometry in the body of the stent within the primary vessel.

6.4 Stent System Removal – Pre-deployment

- If unusual resistance is felt at any time during lesion access before stent implantation, the stent system and the guide catheter should be removed as a single unit (see note below for instructions on this removal method).
- Retraction of an unexpanded stent back into the guide catheter could result in stent or coating damage or stent dislodgment from the balloon. If retraction of the unexpanded stent back into the guide catheter is required, ensure that the guide catheter is coaxially aligned with the stent system and cautiously withdraw the stent system into the guide catheter under direct fluoroscopic visualization.
- Stent retrieval methods (use of additional wires, snares, and/or forceps) may result in additional trauma to the vascular site. Complications can include bleeding, hematoma, or pseudoaneurysm.

Note: When removing the entire stent system and guide catheter as a single unit the following steps should be executed in the order indicated under direct visualization using fluoroscopy.

- If greater than usual resistance is felt during delivery system withdrawal, pay particular attention to guide catheter position. In some cases it may be necessary to pull back slightly on the guide catheter in order to prevent deep seating (unplanned advancement) of the guide catheter and subsequent vessel damage. In cases where unplanned guide catheter movement has occurred, angiographic assessment of the coronary tree should be undertaken to ensure that there is no damage to the coronary vasculature.
- Maintain guidewire placement across the lesion during the entire removal process. Carefully pull back the stent system until the proximal balloon marker of the stent system is just distal to the guide catheter distal tip.
- The stent system and the guide catheter should be pulled back as a single unit until the tip of the guide catheter is just distal to the arterial sheath, allowing the guide catheter to straighten. Carefully retract the un-deployed stent into the tip of the guide catheter and remove the stent system and the guide catheter from the patient again as a single unit while leaving the guidewire across the lesion.

Failure to follow these steps, and/or the use of excessive force to the stent system can potentially result in stent or coating damage, stent dislodgment from the balloon, and/or damage the delivery system.

6.5 Stent System Removal – Post-deployment

- Following stent placement, confirm complete balloon deflation. Ensure balloon is fully deflated before delivery system withdrawal. Larger and longer balloons will take more time to deflate than smaller and shorter balloons. Deflation time is ≤ 30 seconds. Before withdrawing the stent delivery system visually confirm complete balloon deflation under fluoroscopy.
- If greater than usual resistance is felt during delivery system withdrawal, pay particular attention to guide catheter position. In some cases it may be necessary to pull back slightly on the guide catheter in order to prevent deep seating (unplanned advancement) of the guide catheter and subsequent vessel damage. In cases where unplanned guide catheter movement has occurred, angiographic assessment of the coronary tree should be undertaken to ensure that there is no damage to the coronary vasculature.
- If greater than usual resistance is felt during delivery system withdrawal into the guide catheter, see **Stent System Removal - Pre-deployment** for directions.

6.6 Post-Procedure

Care must be exercised when crossing a newly deployed stent with any wire, catheter or ancillary device to avoid disrupting the stent placement, apposition, geometry, and/or coating.

6.7 Brachytherapy

The safety and effectiveness of the Promus ELITE™ Stent in patients with prior brachytherapy of the target lesion have not been established. The safety and effectiveness of the use of brachytherapy to treat in-stent restenosis in the Promus ELITE Stent have not been established. Both vascular brachytherapy and the Promus ELITE Stent alter arterial remodeling. The interaction, if any, between these two treatments has not been determined.

6.8 Magnetic Resonance Imaging (MRI)

Through non-clinical testing, the Promus ELITE Stent has been shown to be MR Conditional (poses no known hazards under specified conditions.) The conditions are as follows:

- Field strengths of 3 Tesla or less with
- Static magnetic field gradient < 14 T/m (extrapolated)
- Product of static magnetic field and static magnetic field gradient < 25 T²/m (extrapolated)
- A calculated rate of change of magnetic field (dB/dt) of 60 T/s or less
- A maximum whole body averaged specific absorption rate (SAR) of lower than 2.0 W/kg for a total active MR scan time (with RF exposure) of 15 minutes or less

The Promus ELITE Stent should not migrate in this MRI environment. MR imaging within these conditions may be performed immediately following the implantation of the stent. This stent has not been evaluated to determine if it is MR Conditional beyond these conditions.

3.0 Tesla Temperature Information

Non-clinical testing of RF-induced heating was performed at 123 MHz in a 3.0 Tesla Magnetom Trio™, Siemens Medical Solutions MR system, software version Numaris/4, syngo™ MR A30A. The stents were in a location and orientation in the phantom that produced the worst case Radio Frequency (RF) heating. RF power was applied for 15 minutes and the measured conductivity of the phantom material was about 0.50 S/m. The phantom average SAR calculated using calorimetry was 2.3 W/kg. The maximal in-vitro temperature rise was calculated as 2.6°C when the local SAR was scaled to 2.0 W/kg for a measured overlapped stent length of 74 mm. Other stent lengths exhibited a lower temperature rise. Predicted in-vivo heating based on these non-clinical tests and computer simulation of the patient exposure to the electromagnetic fields in MRI yielded to the following maximal in vivo rises: for landmarks at the chest level, the calculated temperature rise was 2.6°C with a calculated uncertainty upper bound temperature of 4.7°C for a whole body average SAR value of 2.0 W/kg and a continuous scan time of 15 minutes.

The actual in vivo rise is expected to be less than these values as the calculations did not include the cooling effects due to blood flow in the lumen of the stent and blood perfusion in the tissue outside the stent.

1.5 Tesla Temperature Information

Non-clinical testing of RF-induced heating was performed at 64 MHz in a 1.5 Tesla Intera™ Philips Medical Systems, software version Release 12.6.1.3, 2010-12-02 whole body coil MR scanner. The stents were in a location and orientation in the phantom that produced the worst case RF heating. RF power was applied for 15 minutes and the measured conductivity of the phantom material was about 0.51 S/m. The phantom average SAR calculated using calorimetry was 2.1 W/kg. The maximal in-vitro temperature rise was calculated as 2.6°C when the local SAR was scaled to 2.0 W/kg for a measured single stent length of 39 mm. Other stent lengths exhibited a lower temperature rise. Predicted in-vivo heating based on these non-clinical tests and computer simulation of the patient exposure to the electromagnetic fields in MRI yielded to the following maximal in vivo rises: for landmarks at the chest level, the calculated temperature rise was 2.6°C with an uncertainty upper bound temperature of 4.8°C for a whole body average SAR value of 2.0 W/kg and a continuous scan time of 15 minutes.

The actual in vivo rise is expected to be less than these values as the calculations did not include the cooling effects due to blood flow in the lumen of the stent and blood perfusion in the tissue outside the stent.

In vivo, local SAR depends on MR Field strength and may be different than the estimated whole body averaged SAR, due to body composition, stent position within the imaging field, and scanner used, thereby affecting the actual temperature rise. No tests have been performed on possible nerve or other tissue stimulation possible to be activated by strong gradient magnetic fields and resulting induced voltages.

Image Artifact Information

The calculated image artifact extends approximately 8 mm from the perimeter of the device diameter and 5 mm beyond each end of the length of the stent when scanned in non-clinical testing using a Spin Echo sequence. With a Gradient Echo sequence the calculated image artifact extends 7 mm beyond the perimeter of the diameter and 7 mm beyond each end of the length with both sequences partially shielding the lumen in a 3.0 Tesla Intera (Achieve Upgrade), Philips Medical Solutions, software version Release 2.6.3.5 2009-10-12 MR system with a transmit/receive head coil. This testing was completed using ASTM F2119-07 test method.

6.9 Individualization of Patient Treatment

The device carries an associated risk of acute, sub-acute or late thrombosis, vascular complications, and/or bleeding events. Therefore, patients should be carefully selected, and a P2Y₁₂ inhibitor (e.g. clopidogrel, ticlopidine, prasugrel or ticagrelor) must be prescribed post procedure to reduce the risk of stent thrombosis. Aspirin must be administered concomitantly with P2Y₁₂ inhibitor, and then continued indefinitely to further reduce the risk of thrombosis.

Antiplatelet drugs should be used in combination with (Promus ELITE) drug-eluting stents. Physicians should use the information from the large body of clinical evidence for everolimus-eluting stents, coupled with current literature on drug-eluting stents, current European Society of Cardiology recommendation (or other applicable country guidelines) and the specific needs of the individual patient to determine the specific antiplatelet / anticoagulation regimen to be used for their patients in general practice.

It is very important that the patient be compliant with post-procedural antiplatelet recommendations given by their physician. In selected higher risk patients where the physician determines that the risks outweigh the benefits of continued DAPT, it may be reasonable to interrupt or discontinue therapy after 1 month of DAPT based on the low stent thrombosis rates and no observed increased risk of stent thrombosis demonstrated in the current literature. Patients who require premature discontinuation of antiplatelet therapy should be monitored closely and have their antiplatelet therapy restarted as soon as possible per the discretion of their treating physician.

6.10 Drug Interactions

The following describes known drug interactions for orally administered everolimus at significantly higher doses than are present on the Promus ELITE™ Everolimus-Eluting Coronary Stent. Interactions observed at these higher doses may not be relevant to the Promus ELITE Everolimus-Eluting Coronary Stent System.

When taken orally, everolimus is extensively metabolised by the cytochrome P4503A4 (CYP3A4) in the gut wall and liver and is a substrate for the countertransporter P-glycoprotein. Therefore, absorption and subsequent elimination of everolimus may be influenced by drugs that affect these pathways. Concurrent treatment with strong 3A4 inhibitors and inducers is not recommended unless the benefits outweigh the risk. Inhibitors of P-glycoprotein may decrease the efflux of everolimus from intestinal cells and increase everolimus blood concentrations. In vitro, everolimus was a competitive inhibitor of CYP3A4 and of CYP2D6, potentially increasing the concentrations of drugs eliminated by these enzymes. Thus, caution should be exercised when coadministering everolimus with 3A4 and 2D6 substrates with a narrow therapeutic index. Everolimus has also been shown to reduce the clearance of some prescription medications when it was administered orally along with cyclosporine (CsA).

Please Note: Everolimus when prescribed as an oral medication may interact with the following drugs or substances:

- CYP3A4 isozyme inhibitors (ketoconazole, itraconazole, ritonavir, erythromycin, clarithromycin, fluconazole, calcium channel blockers)
- Inducers of CYP3A4 isozyme (rifampicin, rifabutin, carbamazepin, Phenytoin, phenytoin)
- Antibiotics (Ciprofloxacin, ofloxacin)
- Glucocorticoids
- HMGCoA reductase inhibitors (simvastatin, lovastatin)
- Digoxin
- Cisapride (theoretical potential interaction)
- Sildenafil (Viagra™) (theoretical potential interaction)
- Antihistaminics (Terfenadine, astemizole)
- Grapefruit juice

Formal drug interaction studies have not been performed with the Promus ELITE Everolimus-Eluting Coronary Stent System. Therefore, due consideration should be given to the potential for both systemic and local drug interactions in the vessel wall when deciding to place the Promus ELITE Everolimus-Eluting Coronary Stent in a subject taking a drug with known interaction with everolimus.

6.11 Pregnancy

This product has not been tested in pregnant women or in men intending to father children; effects on the developing foetus have not been studied. While there is no contraindication, the risks and reproductive effects are unknown. It is not recommended that the Promus ELITE Stent System be used in women attempting to conceive, or who are pregnant.

6.12 Use of Multiple Stents

Potential interactions of the Promus ELITE Stent with other drug-eluting or coated stents have not been evaluated in vivo.

Patients should be treated with no more than 2 planned Promus ELITE Stents. Additional stents may be placed if bailout stenting is required. The use of multiple drug-eluting stents will expose the patient to larger amounts of drug and polymer.

When more than one stent is required and results in stent-to-stent contact, stent materials should be of similar composition to avoid the possibility of corrosion due to the presence of dissimilar metals in a conducting medium. If more than one Promus ELITE Stent is needed to cover the lesion, it is suggested that, to avoid the potential for gap restenosis, the stents be

adequately overlapped (with a minimum of 2 mm overlap). Placing multiple stents of different metals in contact with each other may increase the potential for corrosion, though in vitro tests to assess stent-to-stent contact using a platinum chromium alloy stent in combination with a 316L stainless steel alloy stent suggest there is no increased risk of corrosion with this pair.

7.0 DRUG INFORMATION

7.1 Mechanism of Action

The mechanism by which the Promus ELITE Stent inhibits neointimal growth has not been established. At the cellular level, everolimus inhibits growth factor-stimulated cell proliferation. At the molecular level, everolimus forms a complex with the cytoplasmic protein FKBP-12 (FK 506 Binding Protein). This complex binds to and interferes with FRAP (FKBP-12 Rapamycin Associated Protein), also known as mTOR (mammalian Target of Rapamycin), leading to inhibition of cell metabolism, growth and proliferation by arresting the cell cycle at the late G1 stage.

7.2 Drug Interactions

See Section 6.10, Drug Interactions.

7.3 Carcinogenicity, Genotoxicity, and Reproductive Toxicology

A 26-week carcinogenicity study was conducted to evaluate the carcinogenic potential of PROMUS™ (Xience V™) Stents following subcutaneous implantation in transgenic mice. During the course of the study, there were no abnormal clinical observations that suggested a carcinogenic effect of the test group PROMUS (Xience V) Stent. The test group did not demonstrate an increased incidence of neoplastic lesions when compared to the negative control group. However, the positive control and the experimental positive control groups demonstrated notable increases in the incidence of neoplastic lesions compared to either the test or the negative control group. Based on the results of this study, the PROMUS (Xience V) Stent does not appear to be carcinogenic when implanted in transgenic mice for 26 weeks.

In addition, a reproductive toxicity (teratology) study was conducted to demonstrate that implantation of PROMUS (Xience V) Stents in female Sprague-Dawley rats does not affect their fertility or reproductive capability and shows a lack of any reproductive toxicity on their offspring. The PROMUS (Xience V) Stent did not affect the fertility or reproductive capability of female Sprague-Dawley rats. There was no statistical difference between the test article PROMUS (Xience V) Stent and the control system in terms of any of the evaluated parameters. The test article had no effect on litter size and caused no increase of in utero mortality. Additionally, the PROMUS (Xience V) Stent did not cause any reproductive toxicity in the offspring in this study.

8.0 ADVERSE EVENTS

Potential adverse events (in alphabetical order) which may be associated with the implantation of a coronary stent in a native coronary artery include those risks associated with percutaneous transluminal coronary angioplasty as well as additional risks related to the use of the stent as listed below.

- Abrupt closure
- Allergic reaction (including to medications, contrast, stent materials)
- Aneurysm (coronary)
- Angina
- Arrhythmias, including ventricular fibrillation and ventricular tachycardia
- Arteriovenous fistula
- Bleeding
- Cardiac tamponade
- Cardiogenic shock
- Death
- Emboli (including air, tissue, thrombus, plaque or device materials)
- Heart failure
- Hematoma
- Hemorrhage
- Hypotension/hypertension
- Infection, local and/or systemic

- Ischemia, myocardial
- Myocardial infarction
- Pain
- Pericardial effusion
- Pseudoaneurysm, femoral
- Pulmonary edema
- Renal insufficiency or failure
- Respiratory failure
- Restenosis of stented segment
- Shock
- Stent embolisation
- Stent fracture
- Stent migration
- Stent thrombosis and/or vessel occlusion
- Stroke/cerebrovascular accident/transient ischemic attack
- Total occlusion of coronary artery
- Vessel spasm
- Vessel injury (including dissection, perforation, rupture or trauma)

Adverse events associated with daily oral administration of everolimus:

- Abdominal pain
- Acne
- Anemia
- Coagulopathy
- Diarrhea
- Edema
- Hemolysis
- Hypercholesterolemia
- Hyperlipidemia
- Hypertension
- Hypertriglyceridemia
- Hypogonadism male
- Leukopenia
- Liver function test abnormal
- Lymphocele
- Myalgia
- Nausea
- Pain
- Pneumonia
- Pyelonephritis
- Rash
- Renal tubular necrosis
- Sepsis
- Surgical wound complication
- Thrombocytopenia
- Urinary tract infection
- Venous thromboembolism
- Viral, bacterial, and fungal infections
- Vomiting
- Wound infection

There may be other potential adverse events that are unforeseen at this time.

9.0 HOW SUPPLIED

Handling and Storage

Keep dry and protect from light. Store at 25°C (77°F); excursions permitted to 15-30°C (59-86°F). Store product in outer carton.

DO NOT REMOVE FROM FOIL POUCH UNTIL READY FOR USE.

THE FOIL POUCH IS NOT A STERILE BARRIER.

Do not store devices where they are directly exposed to organic solvents or ionizing radiation.

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

The foil pouch contains Argon (Ar) as a storage medium.

10.0 OPERATIONAL INSTRUCTIONS

10.1 Inspection Prior to Use

Check foil pouch for “Use By” date. Do not use the product after the “Use By” date. Carefully inspect the foil pouch and the sterile package before opening. If the integrity of the foil pouch or the sterile package has been compromised prior to the product “Use By” date (e.g., damage of the package), contact your local Boston Scientific representative for return information. Do not use if any defects are noted.

10.2 Materials Required (not included in Stent Delivery System package)

Quantity	Material
1	Appropriate guide catheter (see Table 2.1, Promus ELITE™ Stent System Product Description)
2-3	20 ml (cc) syringe
1000 u/500 cc	Normal heparinized sterile saline
1	≤0.014 in (0.36 mm) guidewire
1	Rotating hemostatic valve
1	Diluted contrast medium 1:1 with normal heparinized sterile saline
1	Inflation Device
1	Torque Device
1	Pre-deployment dilation catheter
1	Three-way stopcock
1	Appropriate arterial sheath

10.3 Preparation

10.3.1 Packaging Removal

Step Action

1. Open the outer box to reveal the foil pouch and carefully inspect the foil pouch for damage.
2. Carefully open the foil pouch by tearing along the tear strip as indicated on the foil pouch to access the sterile barrier package containing the stent delivery system.
3. Carefully inspect the sterile barrier package for damage.
4. Carefully peel open the sterile barrier using aseptic techniques and extract the stent delivery system.
5. Carefully remove the stent delivery system from its protective tubing for preparation of the delivery system. When using a Monorail™ system, do not bend or kink proximal shaft during removal.
6. Remove the product mandrel and stent protector by grasping the catheter just proximal to the stent (at the proximal balloon bond site), and with the other hand, grasp the stent protector and gently remove distally.
7. Examine the device for any damage. If it is suspected that the sterility or performance of the device has been compromised, the device should not be used.

10.3.2 Guidewire Lumen Flush

Step Action

1. Flush stent system guidewire lumen with normal heparinised saline using the flushing needle supplied for the Monorail™ Delivery System at the distal end.
2. Verify that the stent is positioned between the proximal and distal balloon markers. Check for bends, kinks and other damage. Do not use if any defects are noted.

Note: Avoid manipulation of the stent during flushing of the guidewire lumen, as this may disrupt the placement of the stent on the balloon.

Note: Stent contact with any fluid is not recommended as there is a possibility of initiating drug release. However, if it is absolutely necessary to flush the stent with saline, contact time should be limited (1 minute maximum).

10.3.3 Balloon Preparation

Step Action

1. Prepare inflation device/syringe with diluted contrast medium.
2. For treatment of occluded vessels, contrast visualization of the distal vessel to confirm position of guidewire within the lumen is recommended.
3. Attach inflation device/syringe to stopcock; attach to inflation port. Do not bend the hypotube when connecting to inflation device/syringe.
4. With tip down, orient stent system vertically.
5. Open stopcock to stent system; pull negative pressure for 15 seconds; release to neutral for contrast fill.
6. Close stopcock to stent system; purge inflation device/syringe of all air.
7. Repeat steps 5 through 7 until all air is expelled. If bubbles persist, do not use product.
8. If a syringe was used, attach a prepared inflation device to stopcock.
9. Open stopcock to stent system.
10. Leave on neutral.

10.3.4 Delivery Procedure

Step Action

1. Prepare the vascular access site according to standard PTCA practice.
2. Predilate the lesion/vessel with appropriate diameter balloon.
3. Maintain neutral pressure on inflation device attached to stent delivery system.
4. Backload stent delivery system onto proximal portion of guidewire while maintaining guidewire position across target lesion.
5. Fully open rotating hemostatic valve to allow for easy passage of the stent and prevent damage to the stent.
6. Carefully advance the stent delivery system into the hub of the guide catheter. Be sure to keep the proximal shaft straight. Ensure guide catheter stability before advancing the stent delivery system into the coronary artery.

Note: If unusual resistance is felt before the stent exits the guide catheter, do not force passage. Resistance may indicate a problem, and use of excessive force may result in stent damage or stent dislodgment from the balloon. Maintain guidewire placement across the lesion, and remove the stent delivery system and guide catheter as a single unit.

7. Advance the stent delivery system over the guidewire to target lesion under direct fluoroscopic visualization. Utilize the proximal and distal radiopaque balloon markers as a reference point. Fully cover the entire lesion and balloon treated area. The stent should adequately cover healthy vessel proximal and distal to the lesion. If the position of the stent is not optimal, it should be carefully repositioned or removed (see also Precautions, Section 6.4, Stent System Removal – Pre-deployment). The inside edges of the marker bands indicate both the stent edges and balloon shoulders. Expansion of the stent should not be undertaken if the stent is not properly positioned in the target lesion segment of the vessel.

Note: If unusual resistance is felt at any time during lesion access before stent implantation, the stent delivery system and the guide catheter should be removed as a single unit. (See also Precautions, Section 6.4, Stent System Removal – Pre-deployment). Once the stent delivery system has been removed do not re-use.

8. Sufficiently tighten the rotating hemostatic valve. The stent is now ready to be deployed.

10.3.5 Deployment Procedure

Step Action

1. Inflate the delivery system expanding the stent to a minimum pressure of 11 atm (1117 kPa). Higher pressure may be necessary to optimize stent apposition to the arterial wall. Accepted practice generally targets an initial deployment pressure that would achieve a stent inner diameter of about 1.1 times the distal reference vessel diameter (see Tables 11.1 and 11.2). Balloon pressure must not exceed rated burst pressure of 18 atm (1827 kPa) for the 2.25 mm – 2.75 mm diameter stents and 16 atm (1620 kPa) for the 3.00 mm - 4.00 mm diameter stent sizes (see Tables 11.1 and 11.2).
2. Maintain inflation pressure for 15-30 seconds for full expansion of the stent.
3. Deflate balloon by pulling negative pressure on inflation device until balloon is fully deflated. Ensure balloon is fully deflated before delivery system withdrawal. Larger and longer balloons will take more time to deflate than smaller and shorter balloons. Deflation time is ≤30 seconds. Before withdrawing the stent delivery system visually confirm complete balloon deflation under fluoroscopy.
4. Confirm stent position and deployment using standard angiographic techniques. For optimal results, the entire stenosed arterial segment should be covered by the stent. Fluoroscopic visualization during stent expansion should be used in order to properly judge the optimum expanded stent diameter as compared to the proximal and distal coronary artery diameter(s). Optimal expansion requires that the stent be in full contact with the artery wall.
5. If stent sizing/apposition requires optimization, readvance the stent delivery system balloon, or another high-pressure, balloon catheter of the appropriate size, to the stented area using standard angioplasty techniques.
6. Inflate the balloon to the desired pressure while observing under fluoroscopy (refer to product labeling and/or Tables 11.1 and 11.2 for balloon compliance chart). Deflate the balloon. Ensure balloon is fully deflated before delivery system withdrawal. Larger and longer balloons will take more time to deflate than smaller and shorter balloons. Deflation time is ≤30 seconds. Before withdrawing the stent delivery system visually confirm complete balloon deflation under fluoroscopy.
7. If more than one Promus ELITE™ Stent is needed to cover the lesion and balloon treated area, it is suggested that, to avoid the potential for gap restenosis, the stents be adequately overlapped. To ensure that there are no gaps between stents, the balloon marker bands of the second Promus ELITE Stent should be positioned inside of the deployed stent prior to expansion.
8. Reconfirm stent position and angiographic result. Repeat inflations until optimal stent deployment is achieved.

10.3.6 Removal Procedure

Step Action

1. Ensure balloon is fully deflated before delivery system withdrawal. Larger and longer balloons will take more time to deflate than smaller and shorter balloons. Deflation time is ≤30 seconds. Before withdrawing the stent delivery system visually confirm complete balloon deflation under fluoroscopy.
2. Fully open rotating hemostatic valve.
3. While maintaining guidewire position and negative pressure on inflation device, withdraw delivery system.
4. Repeat angiography to assess the stented area. If an adequate expansion has not been obtained, exchange back to the original stent delivery catheter or exchange to another balloon catheter of appropriate balloon diameter to achieve proper stent apposition to the vessel wall.

10.4 Post-Deployment Dilatation of Stented Segments

Precaution: Do not dilate the stent beyond the limits tabulated below.

Nominal Stent Diameter (ID)	Dilatation Limits (ID)*
2.25 mm	2.75 mm
2.50 mm – 2.75 mm	3.50 mm
3.00 mm – 3.50 mm	4.25 mm
4.00 mm	5.75 mm

*Max Stent Inner Diameter

All efforts should be taken to assure that the stent is not under-dilated. If the deployed stent size is still inadequate with respect to vessel diameter, or if full contact with the vessel wall is not achieved, a larger balloon may be used to expand the stent. The stent may be expanded using a low profile and high pressure balloon catheter. If this is required, the stented segment should be recrossed carefully with a prolapsed guidewire to avoid dislodging the stent. The balloon should be centered within the stent and should not extend outside of the stented region. For in-stent restenosis, where details of the original stent are known, the expanded inner diameter of the new stent should not exceed the dilation limits of the original stent. Where details of the original stent are not known, the expanded inner diameter of the new stent should not exceed the reference vessel diameter.

Note: In line with Section 6.6, Post-Procedure: Care must be exercised when crossing a newly deployed stent with any wire, catheter or ancillary device to avoid disrupting the stent placement, apposition, geometry, and/or coating.

11.0 IN VITRO INFORMATION

11.1 Promus ELITE™ Stent System Compliance

Table 11.1. Typical Promus ELITE Stent System Compliance

Pressure atm - kPa	Stent I.D. (mm)						
		2.25	2.50	2.75	3.00	3.50	4.00
8 - 814			2.29	2.50	2.72	3.24	3.72
9 - 910		2.13	2.37	2.58	2.81	3.34	3.81
10 - 1014		2.19	2.43	2.65	2.88	3.43	3.89
11 - 1117	Nominal	2.24	2.50	2.72	2.95	3.51	3.96
12 - 1213		2.29	2.55	2.78	3.01	3.58	4.02
13 - 1317		2.34	2.60	2.84	3.06	3.63	4.08
14 - 1420		2.38	2.65	2.89	3.10	3.68	4.13
15 - 1517		2.42	2.68	2.93	3.14	3.73	4.17
16 - 1620*		2.45	2.72	2.96	3.17	3.77	4.21
17 - 1724		2.47	2.75	2.99	3.20	3.81	4.25
18 - 1827*		2.50	2.77	3.03	3.24	3.85	4.30
19 - 1924		2.52	2.80	3.06	3.28	3.91	4.36
20 - 2027		2.55	2.83	3.09	3.32	3.97	4.43
21 - 2130		2.57	2.87	3.13			
22 - 2227		2.60	2.90	3.17			

*RATED BURST PRESSURE. DO NOT EXCEED.
Nominal Pressure = 11 atm – 1117 kPa

Table 11.2. Typical Promus ELITE Stent System Compliance

Pressure atm - kPa	Stent O.D. (mm)						
		2.25	2.50	2.75	3.00	3.50	4.00
8 - 814			2.48	2.69	2.91	3.43	3.92
9 - 910		2.32	2.56	2.77	3.00	3.53	4.01
10 - 1014		2.38	2.62	2.84	3.07	3.62	4.09
11 - 1117	Nominal	2.43	2.69	2.91	3.14	3.70	4.16
12 - 1213		2.48	2.74	2.97	3.20	3.77	4.22
13 - 1317		2.53	2.79	3.03	3.25	3.82	4.28
14 - 1420		2.57	2.84	3.08	3.29	3.87	4.33
15 - 1517		2.61	2.87	3.12	3.33	3.92	4.37
16 - 1620*		2.64	2.91	3.15	3.36	3.96	4.41
17 - 1724		2.66	2.94	3.18	3.39	4.00	4.45
18 - 1827*		2.69	2.96	3.22	3.43	4.04	4.50
19 - 1924		2.71	2.99	3.25	3.47	4.10	4.56
20 - 2027		2.74	3.02	3.28	3.51	4.16	4.63
21 - 2130		2.76	3.06	3.32			
22 - 2227		2.79	3.09	3.36			

*RATED BURST PRESSURE. DO NOT EXCEED.
Nominal Pressure = 11 atm – 1117 kPa

12.0 CLINICAL TRIALS

The principal safety and effectiveness information for the Promus ELITE Stent System is derived from the global PLATINUM Clinical Trial Program, a series of clinical trials conducted on the PROMUS Element™ Stent System. The PROMUS Element and Promus ELITE Stents utilize the same platinum chromium alloy and the same Everolimus and PVDF-HFP coating, resulting in a similar kinetic release profile. Given the similarities between the PROMUS Element and Promus ELITE Stent Systems and supportive bench and animal study information, the findings from the PLATINUM clinical studies are applicable to the Promus ELITE Stent System.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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Magnetom Trio and syngo are trademarks of Siemens Aktiengesellschaft Corp.
Intera is a trademark of Koninklijke Philips Electronics N.V. Corp.
Viagra is a trademark of Pfizer Inc.
© 2017 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific (Australia) Pty Ltd