# Manuale d'uso Medtrum EasyView

*Per Caregivers*

Humalog® is a trademark of Eli Lilly and Company.

NovoRapid® is a trademark of Novo Nordisk A/S.

Apidra® is a trademark of Sanofi S.A.

Medtrum Technologies Inc.

7F, Building 8, No. 200, Niudun Road

Shanghai 201203, China

Tel: +86-21-50274781

Fax: +86-21-50274779

[www.medtrum.com](http://www.medtrum.com)

Medtrum B.V.

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5683 CJ Best

The Netherlands

Tel: +31 (0) 499745037

**CE 0197**

This product complies with Directive 93/42/EEC (MDD) and Directive 2014/53/EU (RED).

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

The image shows a mobile app sign-up screen titled "Sign up". At the top, there are two radio buttons: "Patient" (unselected) and "Caregiver" (selected). Below this is a dropdown menu showing "United Kingdom (mmol/L)". The form contains four input fields: "Enter email address" (with a "Get code" button to its right), "Enter code", "Enter your full name", and "Enter password". At the bottom is a large blue "Next" button.

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We will only process your Healthcare data through Medtrum App/website with your consent.

If you choose to register your Medtrum account and access our App using the account, the healthcare data recorded on the App will be uploaded to our server in Europe and be held for 2 years. You have the right not to upload any device data to the App, in which case the data will not be accessible to others. If you choose not to register in Medtrum App/website, the only record of your Healthcare data is on the device itself.

If you are under 16 years of age, you will need the consent of your parent or legal guardian to register in this App.

☒ I Agree to [the Terms of Use of Medtrum](#)

☒ I Agree to [the Privacy Policy of Medtrum](#)

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

Connections <span>Add Connection</span>				
Name	SG (mmol/L)	Basal (U/h)	Active Insulin (U)	Action
Alan				Stop following
Voss				Stop following
Eva	8.3		0	Stop following
Wall	6.9	0.8	0	Stop following
Carla	8.1	0.9	0	Stop following
Clair				Stop following

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
Connections		Add Connection		
Name	SG (mmol/L)	Basal (U/H)	Active insulin (U)	Action
 Alan				Stop following
 Voss				Stop following
 Eva	8.3		0	Stop following
 Wall	6.9	0.8	0	Stop following
 Carla	8.1	0.9	0	Stop following
 Clar				Stop following

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
A form titled 'Patient profile' with a light blue header. The form contains the following fields: 'Username:' with the value 'test18@medtrum.cn', 'Photo:' with a person icon, 'Full Name:' with the value 'Clair Simpson', 'Alias:' with an empty text box, 'Mobile:' with a placeholder '--', 'Telephone:' with a placeholder '--', and 'Email:' with a placeholder '--'. At the bottom right, there are two buttons: 'Save' and 'Cancel'.

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
The image shows a 'Report Settings' form for glucose monitoring. At the top, there's a title 'Report Settings (mmol/L)' and 'Save' and 'Cancel' buttons. Below is a 'Show alerts' toggle set to 'On'. The 'Glucose Target Range' section includes 'Overall' (4.5 - 8.8) and 'Hypoglycemia' (3.1). The main section is divided into three columns: 'Glucose Target Range', 'Time Period', and 'Post Meal Analysis'. It lists meal times: Breakfast, Lunch, and Dinner, each with 'Before' and 'After' targets and time ranges. A 'Night-time' row is also present. A 'Restore defaults' button is at the bottom right.

		Glucose Target Range		Time Period		Post Meal Analysis	
		Low	High	From	To	From	To
Breakfast	Before	4.6	6.5	06:15	07:00	01:00	01:30
	After	5.8	8.7				
Lunch	Before	4.6	6.5	11:00	12:00	01:00	01:30
	After	5.8	8.7				
Dinner	Before	4.6	6.5	16:00	17:00	01:00	01:30
	After	5.8	8.7				
Night-time		4.6	6.5	21:00	05:30		

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
BG		CGM		Insulin		Carbs		Exercise	
Average		Average		Average Daily Dose		Average Carbs/Day		Average time/Day	
9.4 mmol/L		9.1 mmol/L		37.33 U		44 g		--	
SD	#	SD	#	SD	# days	SD	#	SD	#
3.17	30	3.46	4901	2.57	4	18.6	6	--	--

Estimated A1C: 7.3% calculated from SG values

Note: Estimated A1C does not replace lab measurement

The second part displays the percentile report, in which you can view the SG distribution at certain daily time points during the selected days.

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
Delete Create		Period: 12/11/2019 - 18/11/2019, All					Set	Print
	Time	Event Type	Details	Notes	Source	Action		
<input type="checkbox"/>	18-11-2019 21:29:00	BG	12.8 mmol/L		PDM	Edit Delete		
<input type="checkbox"/>	18-11-2019 21:29:00	Carbs	66 g		PDM	Edit Delete		
<input type="checkbox"/>	18-11-2019 21:29:00	Exercise	25 min Light		PDM	Edit Delete		
<input type="checkbox"/>	18-11-2019 13:49:00	Injection	10 U		app	Edit Delete		
<input type="checkbox"/>	18-11-2019 13:49:00	Injection	10 U		app	Edit Delete		

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



A screenshot of the 'Quick Print' settings form. The form has a title 'Quick Print' and a subtitle 'Quick Print Settings:'. It includes a 'Select duration:' dropdown set to 'One Week', an 'End date:' dropdown set to 'Yesterday', and a 'Period:' field showing '29/10/2020 - 04/11/2020'. Below this is a 'Week Filter:' section with radio buttons for 'Weekdays' and 'Weekends', and checkboxes for 'Mon', 'Tue', 'Wed', 'Thu', 'Fri', 'Sat', and 'Sun'. The 'Select reports to print:' section has a list of reports with checkboxes: 'All', 'Dashboard', 'Daily Summary', 'Day by Day Overview', 'Sensor Overlay', 'Trend Analysis', and 'Alert'. Each report has a 'From' and 'To' date range with a calendar icon. At the bottom right is a blue button labeled 'Create Report'.

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

A screenshot of the 'My profile' section in the Medtrum app. At the top, there's a header with 'My profile' on the left and a 'Home' button with a right arrow on the right. Below the header, the form contains the following fields: 'Username:' with the value 'test22@medtrum.cn'; 'Profile photo:' with a circular profile picture icon; 'Full Name' with the value 'test22'; 'Gender' with a dropdown menu showing 'Unspecified'; 'Date of birth' with a date picker icon; 'Address' with two stacked text input fields; 'Postcode' with a text input field; 'Mobile' with a text input field; 'Telephone' with a text input field; 'Email' with a text input field; and 'Bio' with a text input field. At the bottom right of the form, there are two buttons: 'Save' and 'Cancel'.

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