# Manuale d'uso per Medtrum EasyView PRO

*Personale sanitario*

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
	Name	Gender	Age	Type	Duration	Account creator	Upload time
<input type="checkbox"/>	Linda	Unspecified	47	Unspecified	47 year 10 month	Crisian	28-08-2018 10:09:44
<input type="checkbox"/>	George	Unspecified	29	Predabetes	47 year 10 month	Crisian	03-08-2018 11:16:10
<input type="checkbox"/>	Green	Unspecified	47	Unspecified	47 year 10 month	Crisian	
<input type="checkbox"/>	Jack	Unspecified	47	Unspecified	47 year 10 month	Crisian	
<input type="checkbox"/>	buhr	Unspecified		Unspecified		Crisian	30-10-2018 11:04:24

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
The 'New Patient' form is a web-based interface for adding a new patient. It features a header with a back arrow and the title 'New Patient', and a 'Home' link in the top right. The form contains several input fields: 'Full Name' (text), 'Gender' (dropdown), 'Height' (text with 'cm' unit), 'Weight' (text with 'kg' unit), 'Date of birth' (date picker), 'Type of diabetes' (dropdown), 'Date of diagnosis' (date picker), 'Address' (two stacked text fields), 'Postcode' (text), 'Mobile' (text), 'Telephone' (text), and 'Email' (text). At the bottom, there are three fields: 'Inpatient' (dropdown), 'Patient no.' (text), and 'Bed no.' (text). The form concludes with 'Create' and 'Cancel' buttons.

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The 'Find Patient' form is a simple search interface. It has a header with a back arrow and the title 'Find Patient', and a 'Home' link in the top right. The main area contains a single text input field with the placeholder text 'Enter username (email)'. Below the input field is a blue 'Search' button.

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
[Create](#) [Add Patient](#) [Find](#)

	Name	Gender	Age	Type	Duration	Account creator	Upload time
<input type="checkbox"/>	 Linda	Unspecified	47	Unspecified	47 year 10 month	Clinician	38-06-2017 10:08:44
<input type="checkbox"/>	 George	Unspecified	26	Predibeter	47 year 10 month	Clinician	03-06-2017 10:30:10
<input type="checkbox"/>	 Green	Unspecified	47	Unspecified	47 year 10 month	Clinician	
<input type="checkbox"/>	 Jack	Unspecified	47	Unspecified	47 year 10 month	Clinician	
<input type="checkbox"/>	 Baker	Unspecified		Unspecified		Clinician	30-10-2017 10:54:24

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

A screenshot of the 'Patient Profile' form. At the top is a navigation bar with icons for Patient List, gao4, Settings, Monitor, Reports, Records, and Quick Print. The main form area is titled 'Patient Profile' and contains the following fields: Username (gao4), Profile photo (a person icon), Full name (gao4), Gender (Unspecified), Height (cm), Weight (kg), Date of birth, Type of diabetes (Unspecified), Date of diagnosis, Address, Postcode, Mobile, Telephone, Email, Inpatient (dropdown), Patient no. (text box), and Bed no. (text box). At the bottom right are 'Save' and 'Cancel' buttons.

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A screenshot of the 'Report Settings (mmol/L)' form. At the top, there's a 'Report Settings' tab and a 'Data to view' dropdown. Below the title bar, there are 'Save' and 'Cancel' buttons. The form is divided into sections. The first section is 'Show alert' with 'On' and 'Off' buttons. The second section is 'Glucose Target Range' with 'Overall' and 'Hypoglycemia' settings. The third section is 'Glucose Target Range' with 'Time Period' and 'Post Meal Analysis' sub-sections. The 'Time Period' section has a table with rows for 'Breakfast', 'Lunch', 'Dinner', and 'Night-time', each with 'Before' and 'After' columns and 'Low' and 'High' target ranges. The 'Post Meal Analysis' section has a table with rows for 'Breakfast', 'Lunch', and 'Dinner', each with 'From' and 'To' columns and 'Hours'.

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Glucose Statistics


BG		CGM		Insulin		Carbs		Exercise	
Average		Average		Average Daily Dose		Average Carbs/Day		Average time/Day	
9.4 mmol/L		9.1 mmol/L		37.33 U		44 g		---	
SD	#	SD	#	SD	# days	SD	#	SD	#
3.17	30	3.46	4901	2.57	4	18.6	6	---	---

Estimated A1C: 7.3% calculated from SG values


Note: Estimated A1C does not replace lab measurement

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
A dialog box titled "Set period" with a close button (X) in the top right corner. It contains a "Duration" dropdown menu set to "7 days". Below this, there are "From" and "To" date fields. The "From" field is set to "19-10-2017" and the "To" field is set to "25-10-2017" with a calendar icon. At the bottom, there are two buttons: "Apply" and "Cancel".

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



A dialog box titled "Set period" with a close button (X) in the top right corner. It contains a "Duration" dropdown menu set to "7 days". Below it, there are "From" and "To" date fields. The "From" field is set to "19-10-2017" and the "To" field is set to "25-10-2017" with a calendar icon. At the bottom, there are two buttons: "Apply" and "Cancel".

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
A dialog box titled "Set period" with a close button (X) in the top right corner. It contains a "Duration" dropdown menu set to "7 days". Below this are "From" and "To" date fields. The "From" field is set to "12-11-2019" and the "To" field is set to "18-11-2019" with a calendar icon. Underneath, there are two rows of radio buttons for selecting days: "Weekdays" (Mon, Tue, Wed, Thur, Fri) and "Weekends" (Sat, Sun). At the bottom are "Apply" and "Cancel" buttons.

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
Bolus

Maximum Bolus: 25U	Bolus Increment: 0.05U	Bolus Calculator: On
Combo/Ext Bolus: On	BG Reminder: Off	Active Insulin Time: 3H

IC Ratios(g/U)

Start	48	--	--	--
Carbs	15	--	--	--

Insulin Sensitivity mmol/L

Start	48	--	--	--
BG	2.8	--	--	--

BG Target

From	48	--	--	--
Range	5.6-5.6	--	--	--

Bolus Reminder (Off)

From	--	--	--
To	--	--	--

CGM

Glucose Alerts: On  
Alert Silence: Off

Glucose Limits				
From	48	--	--	--
Low	4.4	--	--	--
High	13.3	--	--	--

Predictive Alerts (h:mm)			Rate Alerts (mmol/L/min)			Snooze (h:mm)		
High	Off	--	Rise	Off	--	High	Off	--
Low	Off	--	Fall	Off	--	Low	Off	--

Calibration		
Calibration Reminder	Off	--
Calibration Repeat	01:00	
Graph Range	24H	

PDM

Audio Options: Audio+Vibrate	Patch Expiration: On 2H	Sensor Expiration: Off
Low Reservoir: Insulin 10U	Auto Off: Off	Confidence Reminder: On
Hourly Max: 25U	Daily Max: 80U	Low Suspend: On 3.9mmol/L
Predictive Low Suspend: On 00:35 min	Out Of Range: --	

Time Format: 24H	Username: 0703664461	Display Settings (Screen): 60sec
Display Settings (Backlight): 60sec	Device Options (CGM System): On	Device Options (Insulin Pump): On
Glucose Unit: mmol/L	Alarm Clock: Off	Telephone: --

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Complications:
Create

Concomitant diseases	Date of diagnosis	Suggestions	Action
----------------------	-------------------	-------------	--------

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

A screenshot of the 'Quick Print' settings interface. At the top, it says 'Quick Print' and 'Quick Print Settings:'. Below this, there's a 'Select duration:' dropdown set to 'One Week'. The 'End date:' is 'Yesterday' and the 'Period:' is '29/10/2020 - 04/11/2020'. Under 'Week Filter:', there are radio buttons for 'Weekdays' (selected), 'Weekends', and 'All days' (with sub-options for Mon, Tue, Wed, Thu, Fri, Sat, Sun). The 'Select reports to print:' section has a list of report types with radio buttons: 'All' (selected), 'Dashboard', 'Daily Summary', 'Day by Day Overview', 'Sensor Overlay', 'Trend Analysis', 'Logbook', 'Event Summary', 'Alert', and 'Device Settings Review'. To the right of each report type is a 'From' and 'To' date range selector, all showing '29-10-2020' to '04-11-2020'. A 'Date' dropdown is next to the 'Event Summary' row. At the bottom right, there is a blue button labeled 'Create Report'.

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


The screenshot shows the 'My profile' page of the Medtrum application. At the top, there is a header with the Medtrum logo on the left and a user profile icon with the name 'Kevin' on the right. Below the header, the page title 'My profile' is followed by a 'Home' link. The main content area contains a form for editing the user's profile. The form fields are: Username (pre-filled with 'Kevin'), Profile photo (with a placeholder icon), Full Name (pre-filled with 'kevin1'), Gender (a dropdown menu currently showing 'Unspecified'), Date of birth (a date picker), Profession (a dropdown menu currently showing 'Unspecified'), Organization, Address (a multi-line text area), Postcode, Mobile, Telephone, Email, and Bio. At the bottom right of the form, there are 'Save' and 'Cancel' buttons.

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The screenshot shows the 'Change Password' page of the Medtrum application. At the top, there is a header with the text 'Change Password' on the left and a 'Home' link on the right. The main content area contains a form for changing the password. The form fields are: Old password, New password, and Re-enter password. Below the Re-enter password field, there is a 'Forgot your password?' link. At the bottom right of the form, there are 'Save' and 'Cancel' buttons.

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Humalog® is a trademark of Eli Lilly and Company.
NovoRapid® is a trademark of Novo Nordisk A/S.
Apidra® is a trademark of Sanofi S.A.