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Produktbeschreibung

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Instructions for use

Product description

Atrauman dressings are non-medicated impregnated dressings for atraumatic wound treatment.

Composition

Support fabric made of hydrophobic polyester tulle, impregnated with a non-medicated ointment based on triglycerides (neutral fats).

The ointment contains (INCI-assigned names):

Caprylic/Capric/Myristic/Stearic Triglyceride;

Bis-Diglyceryl Polyacyladipate-2

Properties and mode of action

The soft, thin support material made of polyester ensures close contact with the whole surface of the wound which is necessary for the drainage of wound exudate. The water-repellent polyester filaments and the smooth surface structure of the tulle counteract adhesion to the wound. The dressing can be removed painlessly from the wound without causing any new damage. Atrauman's non-adhering properties are further enhanced by the ointment base with which the product is impregnated. It has no occlusive effects and does not cause maceration of the underlying skin. Excessive amounts of exudate are drained freely into an absorbent secondary dressing placed on top of the Atrauman dressing. By keeping the edges of the wound soft and supple, the non-medicated, non-sensitising ointment base enhances epithelialisation. The base is non-medicated and contains no paraffins. Because it leaves no residues on removal, it facilitates both surgical procedures and wound cleansing.

Indications

For the treatment of wounds, especially for exuding wounds (abrasions, lacerations, gashes, burns, scalds, acid burns, for covering of donor and recipient sites in skin grafting, in plastic and cosmetic surgery, for nail removal, circumcisions etc). Due to the neutral properties of the non-medicated ointment base, Atrauman is especially useful in dermatology and for patients with sensitive skin, or those sensitive to certain medicaments.

The physician can at his own discretion apply the dressing in combination with topical medication of his choice.

Mode of application

Remove the Atrauman dressing from the sterile peelpack with both cover papers in place and cut it to the size of the wound if necessary. After removal of one of the cover papers, place this side of the dressing onto the wound and remove the second cover. Place a sterile absorbent dressing pad over the Atrauman dressing to absorb exudate. Unless otherwise ordered by the doctor, apply a new Atrauman dressing at each dressing removal.

Note: If there are therapeutic reasons for not changing Atrauman, the ointment base may be absorbed when left in place for a long period. In this case any wound adhesion which may occur can be eliminated by placing a second Atrauman dressing over the first one.

Contra-indications

Do not use Atrauman on any patients who may be allergic to any of its ingredients.

General instructions

Contents sterile unless peel packs are damaged. Store in a horizontal position at room temperature (not above 40 °C).

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