DISTRIBUTED BY:

B | BRAUN
SHARING EXPERTISE

B. Braun Surgical, S.A.

Carretera de Terrassa, 121
08191 Rubí, Spain

Phone +34 93 586 62 00
Fax +34 93 588 10 96

www.bbraun.com

MANUFACTURED BY:

Synthesia, a.s.
Semtin 103,
530 02 Pardubice
Czech Republic

Phone +420 466 821 111
www.synthesia.eu

CE 1023

Cellistyp[®]

Oxidized Resorbable Cellulose

BG	Резорбируем хемостат с оксидирана целулоза	2
CN	可吸收式氧化纤维素止血产品	6
CZ	Vstřebatelné hemostatikum z oxidované celulózy	10
DE	Oxidiertes resorbierbares Cellulose Hämostyptika	14
DK	Oxideret absorberbar cellulose karklemme	18
ES	Hemostático absorbible de celulosa oxidada	22
FI	Häpetettu resorboiva selluloosahemostaatti	26
FR	Hémostatique résorbable de cellulose oxydée	30
GB	US Oxidized resorbable cellulose hemostat	34
GR	Αιμοστατικό οξειδωμένης απορροφήσιμης κυτταρίνης	38
HU	Elszívódó oxidált cellulóz vérzéscsillapító	42
IT	Cellulosa ossidata emostatica riassorbibile	46
NL	Resorbeerbaar hemostatum uit geoxideerde cellulose	50
NO	Resorberbart oksidert cellulose-hemostatikum	54
PL	Wchłaniany hemostatyk zawierający utlenioną celulozę	58
PT	Hemostático reabsorvível de celulose oxidada	62
RO	Hemostatic absorabil din celuloză oxidată	66
SE	Oxiderat återabsorberande cellulosahemostatikum	70
SK	Vstřebateľné hemostatikum z oxidovanej celulózy	74
TR	Emilebilir Oksitlenmiş Selüloz Kanama Durdurucu	78
RU	Рассасывающийся гемостатический материал из окисленной целлюлозы	82

Cellistyp[®]

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18640-033-3/0

Cellistyp®

Topical Absorbable Hemostats, Sterile
FOR SURGICAL USE

Description:

Cellistyp® is a sterile absorbable knitted fabric (Cellistyp® and Cellistyp® D-K) or a sterile absorbable soft layered material (Cellistyp® F), prepared through the controlled oxidation of cellulose. The product's color is white and it can be cut or sutured without fraying.

The product line includes regular density Cellistyp® products, used to control capillary, minor venous and minor arteriolar bleeding. High density Cellistyp® D-K products providing greater thickness and improved endurance, used to control higher volume capillary, venous or arteriolar bleeding and a fibrillar version of the Cellistyp® F product with reduced weight and extremely high flexibility, used to control hemostasis over a large area, surface applications, irregularly shaped bleeding sites or in areas that are difficult to access, as the desired amount can be grasped with forceps and easily placed onto the bleeding site. Any amount of Cellistyp® F can be easily separated, yet the fibers continue to cohere to one another without uncontrolled dispersal over the operative site.

Mode of Action:

When saturated with blood, Cellistyp® swells into a gelatinous mass, still maintaining its original structure. The product then helps in the formation of a clot by initial denaturation of blood proteins, thus supporting local hemostasis and hemorrhage control.

When used properly in minimal amounts, Cellistyp® is absorbed with practically no tissue reaction, depending on the degree of saturation with blood, and the tissue bed, practically without residues. In addition to topical hemostasis, the bacteriostatic and bactericidal properties of Cellistyp® inhibit growth and multiplication of gram positive and gram negative organisms including aerobic and

anaerobic bacteria. Nevertheless, Cellistyp® is not intended to be a substitute for systemically administered therapeutic or prophylactic anti-microbial agents to control or prevent post-operative infections.

Indications:

Cellistyp® is used adjunctively during surgical and minimal invasive procedures to control capillary, minor venous and minor arteriolar bleeding, when conventional hemostatic methods such as ligation are not effective or practical. Cellistyp® can be cut to size for use in endoscopic procedures.

Device Form:

Cellistyp® (regular density knitted fabric) is supplied sterile in double pouch packaging. Cellistyp® D-K (high density knitted fabric) is supplied sterile in double pouch packaging. Cellistyp® F (regular density fibrillar version) is supplied sterile in double pouch packaging.

Instructions for Use:

Sterile techniques should be observed when removing Cellistyp® from the packaging. The outer pouch should be peeled open to allow the inner pouch to drop on a sterile tray or towel. The inner pouch should be peeled open by a scrubbed, gloved and gowned member of the operating team. The product can then be removed using forceps. The required amount of Cellistyp®, adjusted as necessary, is laid on the bleeding site or held firmly against the tissue until bleeding stops. Any unused product should be discarded.

Contraindications:

Topical Absorbable Hemostats should not be used for implantation in bone lesions, such as fractures, because of potential interference with callus formation and a theoretical chance of cyst formation. Topical Absorbable Hemostats should not be used in conjunction with methyl methacrylate adhesives, for example in orthopedic surgery, because their presence may reduce the adhesive strength of the bonding

agent to bone.

Topical Absorbable Hemostats should not be used around the optic nerve and chiasma opticum and should not remain in contact with other neural pathways, where extra pressure induced by the swelling of the absorbable hemostatic material may affect neural pathway function.

Topical Absorbable Hemostats should not be used to control bleeding from large arteries.

Topical Absorbable Hemostats are absorbable hemostats and should not be used for adhesion prevention.

Topical Absorbable Hemostats should not be used on non-hemorrhagic serous oozing surfaces, since body fluids other than whole blood, such as serum, do not react with the product to produce a satisfactory haemostatic effect.

Topical Absorbable Hemostats are not intended to be used as a substitute for careful surgery and the proper use of sutures and ligatures.

Enclosing the product in a contaminated undrained wound may lead to complications and should be avoided.

Topical Absorbable Hemostats should not be used inside the blood circulatory system to avoid the passage of loose fragments of the product through blood salvage system filters. The presence of product fragments may lead to partial or full obstruction of the vasculature system.

When used during vascular surgery, extra care must be taken to avoid penetration of the product into the vascular system. Otherwise, loose fragments of the product may lead to partial or full obstruction of the circulatory system.

Topical Absorbable Hemostats are not designed, sold or intended for any use other than as indicated. It is the physician's responsibility to determine whether any physical impairment of the patient could contraindicate the use of this device.

Warnings:

Topical Absorbable Hemostats should be used by trained and skilled personnel only.

The device is intended for single use only. Do not reuse or re-sterilize. Discard any unused material. Do not use if the gauze is yellow/brown.

Browning of the gauze and a tendency to crumble may indicate damage to the product caused by incorrect storage conditions.

Do not use the device if the package or contents appear damaged in any way.

Wadding or tight packing of Topical Absorbable Hemostats should be avoided, especially within the bony closure of the central nervous system or within other rigid cavities, where extra pressure induced by the swelling mass of the absorbable hemostatic material may interfere with normal function and may even cause tissue death.

If used as a wrap during vascular surgery, do not apply too tightly.

When Topical Absorbable Hemostats are used to achieve hemostasis within, around, or in the proximity of foramina in bones, in areas of bony confine, the spinal cord, or in the proximity of the optic nerve and chiasma opticum, they must always be removed after hemostasis has been achieved. Otherwise, the material will swell and could exert unwanted pressure.

Topical Absorbable Hemostats should be applied dry and should not be moistened prior to application. Make sure to remove as much blood as possible from the bleeding site before applying the product.

Topical Absorbable Hemostats should not be impregnated with anti-infective agents or with other materials such as buffering or hemostatic substances. The hemostatic effect is not enhanced by the addition of thrombin - thrombin activity is destroyed by the low pH of the product.

In general, use a Redon drain in order to present plugging and to ensure safe drainage. The drain must not be placed on or close to the applied Topical Absorbable Hemostats.

Precautions:

Use only as much of the product as necessary for hemostasis, holding it firmly in place until bleeding

stops.

When used in urological cases, extra care must be exercised to prevent the device plugging the urethra, ureter or a catheter.

Topical Absorbable Hemostats should not be used as a surface dressing except for immediate bleeding control as they inhibit the growth of new skin.

If Topical Absorbable Hemostats are used temporarily to line the cavity of large open wounds, they should be placed so as not to overlap the skin edges. Topical Absorbable Hemostats should not be used in infected areas.

The use of Topical Absorbable Hemostats should not be preceded by the application of escharotic chemicals, such as silver nitrate or other similar products, since absorption of the device may be reduced in chemically cauterized areas.

Potential Complications:

Potential complications include, but are not limited to:

- Stenosis
- Paralysis, nerve damage and / or blindness
- Drainage complications after cholecystectomy and urethral procedures
- "Burning" and "stinging" sensations, sneezing
- Risks normally associated with interventional procedures.

There have been reports of a stenotic effect when oxidized cellulose has been applied as a wrap during vascular surgery. Although it has not been established that the stenosis was directly related to the use of oxidized cellulose, it is important to be cautious and to avoid applying the material tightly as a wrapping.

Paralysis and nerve damage have been reported when oxidized cellulose was used around, within, or in the proximity of foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasma opticum. While most of these reports have been in connection with laminectomy, reports of paralysis have also been received in connection with other procedures. Blindness has been reported

in connection with surgical repair of a lacerated left frontal lobe when oxidized cellulose was placed in the anterior cranial fossa.

Possible prolongation of drainage in cholecystectomy and difficulty in passing urine via the urethra after prostatectomy have been reported. There has been one report of a blocked ureter after kidney resection, in which postoperative catheterization was required.

Occasional reports of "burning" and "stinging" sensations and sneezing when oxidized cellulose has been used as packing in epistaxis are believed to be due to the low pH of the product.

A burning sensation has been reported when oxidized cellulose was applied after nasal polyp removal and after hemorrhoidectomy. Headache, burning, stinging and sneezing in epistaxis and other rhinology procedures, as well as stinging where oxidized cellulose was applied on surface wounds (varicose ulcerations, derm abrasions and donor sites) have also been reported.

Interactions:

Not applicable.

Overdose/intoxication:

Not applicable.

Effects on the ability to drive and operate machinery:

Not applicable.

Use in pregnancy and lactation:

Not applicable.

Storage:

Keep dry.

For storage temperature conditions please refer to the label on the box.

Protect from direct sunlight. Keep in the original packaging.

Do not use beyond the use by date stated on the label. Browning of the gauze and a tendency to

crumble indicate incorrect storage.

Sterilized by gamma irradiation. Do not resterilize.

Any adverse effects not noted here should be reported.

REF CATALOGUE NUMBER

LOT BATCH CODE

SIZE PRODUCT SIZE

STERILE R STERILIZED USING IRRADIATION

USE BY DATE

CAUTION, CONSULT ACCOMPANYING DOCUMENTS

DO NOT RESTERILIZE

DO NOT REUSE

TEMPERATURE LIMITATION (STORAGE REQUIREMENTS)

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